

# Venous Wounds –

## Continuing Clinical Challenges





Up to 6 million people in the United States suffer from venous leg ulcers annually<sup>1</sup> including at least 2.5 percent of the long term care population.<sup>2</sup> In the UK, district nurses spend about 50 percent of their time treating leg ulcers.<sup>3</sup>

# Venous Insufficiency Leads to Venous Wounds

Patients are at risk for venous ulcers and venous dermatitis due to inflammation, hypoxia and edema caused by underlying venous hypertension (abnormally high venous pressure)<sup>4,5</sup> Venous hypertension is caused by chronic venous insufficiency (CVI)<sup>4,5</sup> CVI is the chronic inability to pump enough blood from the legs back to the heart.

CVI is the result of one or more of these pathologies:<sup>6</sup>

- Obstruction of the veins (usually thrombosis, but also increased abdominal pressure from obesity, pregnancy, etc.)
- Failure of the muscle pump related to inactivity, paralysis, decreased ankle range of motion or disease
- Incompetent one-way valves in the veins



**Signs and symptoms associated with late stage Venous Hypertension or CVI include:**

- Venous leg ulcers:<sup>7,8</sup>
  - painful in 90 percent of all cases
  - recurrent in up to 72 percent of all patients
  - at least 50 percent take more than a year to close, with about 34 percent taking more than 5 years
- Lower leg dermatitis:<sup>7,8</sup>
  - itching, with scratching, leading to ulcer formation
  - dry scaling or weeping crusts
  - brown staining (hemosiderin deposits)
  - thickened skin with scales
  - fibrotic skin that causes the lower 1/3 of the lower leg to be thin, giving the appearance of an inverted champagne bottle (lipodermatosclerosis)
  - hypersensitivity to allergens, such as products initially used to decrease itching and dryness or to treat ulcers and infection

## **PolyMem<sup>®</sup>** The ideal venous wound management tool!<sup>™</sup>

**Built right into each PolyMem dressing, these four capabilities are ready when you need them – without incurring extra costs or gathering additional supplies.**

**CLEANSSES** Contains a mild nonionic, nontoxic, tissue friendly cleansing agent, activated by moisture, that is gradually released into the wound bed to promote effective autolytic debridement. The built-in continuous cleansing capabilities usually eliminate the need to cleanse the venous ulcer so you can avoid disrupting the growth of healthy new tissue, cooling the wound by rinsing or causing pain during the dressing change process.

**FILLS** Gently expands to fill and conform to the contours of the venous ulcer and the patient's body, which helps maintain a moist wound healing environment.

**ABSORBS** In order to accommodate the full range of exudate levels, PolyMem wicks up to ten times its weight in exudate from venous ulcers.

**MOISTENS** Keeps the wound bed moist and soothes traumatized tissues, helping to relieve wound pain and providing comfort at the often painful venous ulcer site. The moisturizer also helps keep the dressing pad from adhering to the wound so it removes with virtually no pain or trauma, improving caregiver-patient interaction and the overall care experience.

# PolyMem® The perfect fit for your venous insufficiency patients.

## Why PolyMem for Venous Ulcers and Venous Dermatitis?

- PolyMem helps promote healing of venous ulcers and venous dermatitis.<sup>10,11</sup>
- PolyMem contains a mild nonionic, nontoxic, tissue friendly cleansing agent, which helps to continuously remove the fibrous slough that is a hallmark of venous ulcers. This built in cleanser usually eliminates the need for painful and disruptive manual wound bed cleansing during dressing changes.<sup>10,11</sup>
- PolyMem helps decrease inflammation.<sup>12</sup> Prolonged local inflammation is a key underlying cause of both venous ulcers and venous dermatitis.<sup>3,4,5</sup> Inflammation is the single greatest cause of wound pain.<sup>13,14</sup>
- PolyMem contains a superabsorbent starch copolymer to lock moisture in the dressing and remove exudate from the wound. This can eliminate maceration.<sup>10,11,15,16,17,22,23</sup> PolyMem is available in multiple configurations which can be used to manage venous wounds that range from dry to heavily exuding.
- PolyMem helps relieve the often chronic and severe pain of venous ulcers.<sup>10,15,16,17</sup> PolyMem helps reduce wound pain by altering the actions of certain pain sensing nerve endings.<sup>12</sup>
- PolyMem contains glycerol which helps assure the dressing is non-adherent to the wound bed<sup>18</sup> and reduces the foul odors often associated with venous ulcers.<sup>19</sup>
- PolyMem helps relieve pruritis (itching),<sup>10</sup> which is common with venous dermatitis.<sup>7</sup>
- PolyMem may be used under compression and does not have a tendency to create pressure point injury due to its all around soft structure and conformability to the patient's body.<sup>21</sup>
- PolyMem Silver® contains small particle silver that helps kill absorbed microorganisms without leaching into the wound bed.<sup>20</sup>

## PolyMem's formulation improves venous wound outcomes.

### EXAMPLE OF ACTUAL CLINICAL RESULT - NON-HEALING VENOUS ULCER<sup>23</sup>



Granulation at first weekly dressing change of this highly exuding, venous ulcer; compressed with multi-layered wrap



Second weekly dressing change; decreased edema, closing quickly; multi-layered compression continued



Fourth and final visit, closed. Fitted for compression hosiery

### EXAMPLE OF ACTUAL CLINICAL RESULT - VENOUS DERMATITIS<sup>10</sup>



Fibrin, slough-filled venous dermatitis wounds. Pain, initially 10 (on 0-10 scale), almost completely eliminated by day 5; compressed with multi-layered wrap



Almost closed at six weeks. Pruritis (itching) significantly diminished; multi-layered compression continued




Six months after patient's final wound care visit, skin remains intact


**PolyMem is optimized for each compression type and exudate level.**



		Compression Types				
		Elastic Compression		Inelastic Compression		No Compression
		Multi-layer Wraps	Compression Hosiery	Unna's Boot	Short Stretch Bandages	Stretch Gauze, Etc.
<b>Wound Condition</b>	<b>Heavily exudating venous ulcers</b>	PolyMem Max, <sup>®</sup> or PolyMem Wic, <sup>®</sup> which allows exudate to wick away from skin into absorptive layers		PolyMem Max; add PolyMem Wic under PolyMem Max for extra absorbency		
	<b>Low or moderately exudating venous ulcers</b>	PolyMem; or PolyMem Max when longer wear times are desired				
	<b>Venous dermatitis</b>	PolyMem				
Configurations that include silver are ideal when antimicrobial activity is desired.						



PolyMem and PolyMem Silver



PolyMem Max and PolyMem Silver Max



PolyMem Wic and PolyMem Silver Wic

- Leach MJ. Making sense of the venous leg ulcer debate, a literature review. *Journal of Wound Care* 2004 Feb;13(2): 52-6.
- Wipke-Tevis DD, Rantz MJ, Mehr DR, Popejoy L, Petroski G, Madsen R, Conn VS, Grando VT, Porter R, Maas M. Prevalence, incidence, management and predictors of venous ulcers in the long-term-care population using the MDS. *Advances in Skin and Wound Care* 2000 Sep-Oct;13(5):218-224.
- Jankunas V, Bagdonas R, Samsanavicius D, Rimdeika R. An analysis of the effectiveness of skin grafting to treat chronic venous leg ulcers. *Wounds* 2007 May;19(5):128-137.
- Bergan JJ, Schmid-Schonbein GW, Coleridge Smith PD, Nicolaidis AN, Boisseau MR, Eklof B. Chronic Venous Disease. *New England Journal of Medicine* 2006 Aug;355:488-98.
- Chen WY, Rogers AA. Recent insights into the causes of chronic leg ulceration. *Wound Repair and Regeneration* 2007 Jul-Aug;15(4):434-49.
- Orstead HL, Radke L, Gorst R. The impact of musculoskeletal changes on the dynamics of the calf muscle pump. *Ostomy Wound Management* 2001 Oct;47(10):18-24.
- Valencia IC, Falabella A, Kirsner RS, Eaglstein WH. Chronic Venous Insufficiency and venous leg ulceration. *Journal of the American Academy of Dermatology* 2001 March;44(3):401-21.
- Hofman D, Ryan TJ, Arnold F, Cherry GW, Lindholm C, Bjellerup M, Glynn C. Pain in venous leg ulcers. *Journal of Wound Care* 1997 May;6(5):222-4.
- Fluger SL, Clark RA. Stasis Dermatitis. *Emedicine from WebMD*. <http://www.emedicine.com/derm/topic403.htm>. Reviewed November 29, 2007.
- Harrison J. Persistent multiple venous hypertension dermatitis ulcerations closed completely using polymeric membrane wound filler under compression. Presented at WOCN Society 39th Annual Conference. Poster #1228, June 10 - 13, 2007. Salt Lake City, UT USA.
- Caras J. Three malleolus wounds of two years' duration closed in four weeks using silver polymeric membrane dressings. Presented at 19th Annual Symposium on Advanced Wound Care (SAWC). Poster #82, Apr 30 - May 3, 2006. San Antonio, TX USA.
- Beltz AJ, Newman A, Kahn AR, Ruggles T, Ekmejer L. A Polymeric Membrane Dressing With Antinociceptive Properties: Analysis With a Rodent Model of Stab Wound Secondary Hyperalgesia. *The Journal of Pain*, February, 2004; 5(1):38-47.
- Abraham SE. Pain Management in wound care. *Podiatry Management*. June/July 2006:165-168.
- Levine JD, Reichling DB. Chapter 2 Peripheral Mechanisms of Inflammatory Pain. In Wall PD, Melzak R, Editors. *Textbook of Pain*. 4th edition. Edinburgh, UK: Churchill Livingstone, 1999; pages 59-84.
- Thurs K. Chronic leg ulcer healed in only seven dressing changes using silver polymeric membrane dressings. Presented at Wound Ostomy Continence Nursing (WOCN) Society 39th Annual Conference. Poster #1240, June 10 - 13, 2007. Salt Lake City, UT USA.
- Harrison J. Chronic venous ulcer closing steadily with complete elimination of wound pain using standard or silver polymeric membrane wound filler under compression. Presented at WOCN Society 39th Annual Conference. Poster #1229, June 10 - 13, 2007. Salt Lake City, UT USA.
- Hubbard M. Pain relief and healing using Polymeric Membrane Dressings under compression for venous hypertension ulcers. Presented at 20th Annual Symposium on Advanced Wound Care (SAWC). Poster #160, April 28-May 1, 2007. Tampa, FL USA.
- Foresman PA, Etheridge CA, Rodeheaver G. A Wound Dressing Evaluation on Partial-Thickness Rat Wounds. Symposium on Advanced Wound Care Health Management Publications, Inc., 1991 Annual Meeting. Poster Presentation.
- Vandeputte J. Clinical Relevant Discoveries beyond occlusion when using a glycerin hydrogel. Poster Presented at 1995 Advanced Wound Care Symposium.
- Burd A, Kwok CH, Hung SC, Chan HS, Gu H, Lam WK, Huang L. A comparative study of the cytotoxicity of silver-based dressings in monolayer cell, tissue explant, and animal models. *Wound Repair & Regeneration* 2007 Jan-Feb; 15(1):94-104.
- Reaves K. Only two silver polymeric membrane dressings needed to close two venous ulcers in a diabetic patient. Presented at 22nd Annual Clinical Symposium on Advances in Skin & Wound Care (CSASWC). Poster #155, Oct 11 - 14, 2007. Nashville, TN USA.
- Fowler E, Papan JC. Clinical evaluation of a polymeric membrane dressing in the treatment of dermal ulcers. *Ostomy Wound Manage* 1991; 35:35-8, 40-4.
- Harrison J. Polymeric membrane wound filler under compression: Venous hypertension ulcer closed with only four dressing changes. Presented at Wound Ostomy Continence Nursing (WOCN) Society 39th Annual Conference. Poster #1230, June 10 - 13, 2007. Salt Lake City, UT USA.

**PolyMem<sup>®</sup>**

**Ferris Mfg. Corp.**

5133 Northeast Parkway  
Fort Worth, TX 76106-1822 U.S.A.

Toll Free U.S.A.: 800-765-9636

International: +1 630-887-9797

[www.PolyMemShapes.com](http://www.PolyMemShapes.com)

[www.PolyMem.com](http://www.PolyMem.com)

[www.PolyMem.eu](http://www.PolyMem.eu)

This document is meant for general informational purposes only. See individual product literature for specific indications and instructions for use.

PolyMem, PolyMem Silver, PolyMem Wic, PolyMem Wic Silver, PolyMem Max, PolyMem Max Silver, Shapes, Shapes by PolyMem, The Shape of Healing, Ferris and FMC Ferris and design are trademarks of Ferris. The marks may be registered or pending in the US Patent and Trademark Office and in other countries. © 2012 Ferris Mfg. Corp. All rights reserved. MKL-281,REV-5,0912