




CASE STUDY

Huge Sacral Pressure Ulcer Closed in Four Months Using PolyMem Silver[®] and PolyMem Wic Silver[®] Cavity Filler Dressings



BEFORE



4 MONTHS AFTER TREATMENT

CASE STUDY

Huge Sacral Pressure Ulcer Closed in Four Months Using PolyMem Silver and PolyMem Wic Silver Cavity Filler Dressings

Dr. Charalambos Agathangelou, MD, PhD Geriatrics, Dhali Community Geriatric Home, CYPRUS



SEPT 1, 2006

Initial Assessment in Hospital. Necrotic hard tissue extended deep into the muscle. Wound deteriorated using IV antibiotics & povidone ointment. Pain: 10.



SEPT 13, 2006

After surgical debridement. 4+ pseudomonas. Began PolyMem Wic Silver cavity filler changed twice daily. No systemic medications. Pain: 7 - 8.



SEPT 30, 2006

After only 2½ weeks, yellow slough and exudate are already reduced without additional cleansing or rinsing. Dressing changes are now daily. Pain is 0!

PROBLEM

An alert 85-year-old diabetic woman was diagnosed with a pressure ulcer five days after she was admitted to the hospital for pneumonia. Her wound pain was 10 on the faces scale despite analgesics. The nurses at the hospital were using povidone iodine 10% ointment on the ulcer and were unwilling to change to another treatment. The wound deteriorated; when the patient was dismissed to the geriatric nursing home one-half of her gluteus maximus muscle was necrotic and had to be surgically removed. Her pain was now 7 – 8 on the faces scale.

The geriatric nursing home placed the patient on an airbed mattress. Her gaping yellow-slough-lined strong smelling wound with 4+ pseudomonas needed quick healing to minimize the patient's discomfort as well as treatment costs. Considering her comorbidities, avoiding the use of systemic antibiotics and their inherent side effects was important as well.

RATIONALE

The author chose PolyMem dressings and cavity filler because they are not ordinary foam – they have several key integral components which draw and concentrate healing substances from the body into the wound bed to promote rapid healing. The built-in wound cleanser facilitates autolytic debridement directly by loosening the bonds between the slough and the wound bed. No manual wound cleansing is usually needed, allowing for less disruption of the new growth at the wound bed and very quick and easy dressing changes. The glycerol in the dressing prevents sticking and protects the periwound from maceration. Starch copolymers give the dressings superior absorption, and a thin semi-permeable film backing optimizes moisture and protects the wound. Together these components also allow PolyMem dressings to hydrate the wound bed when needed and help decrease pain. PolyMem dressings also insulate the wound, helping maintain a steady warm temperature, which speeds healing.

PolyMem Silver dressings are effective against pseudomonas, so we felt using this product could keep the patient from requiring antibiotics. Recently several other modern silver dressings were shown to be severely cytotoxic in vivo, but cells in contact with PolyMem Silver dressings proliferated. This further affirms the author's decision to use PolyMem Silver dressings on this patient's wound.



OCT 14, 2006

Cultures are negative after one month of using only PolyMem dressings. The wound is much smaller, as well, with granulation instead of slough.



NOV 11, 2006

Wound continues to decrease in depth, undermining and circumference using PolyMem Silver dressings changed daily with no cleansing.


METHODOLOGY

Layers of PolyMem Wic Silver cavity filler were placed into the wound bed and covered with a PolyMem Silver secondary dressing to prevent leakage. These were changed twice daily at first, then daily, then every three days. No systemic or topical antibiotic therapy was added. For one month after the wound closed, PolyMem dressings were continued on the scarred area to reduce pain and protect the fragile new soft tissue. These dressings were changed every 4 – 5 days.



NOV 28, 2006

Wound is now shallow – PolyMem Wic Silver cavity filler is discontinued. PolyMem Silver dressing is now changed every two days.



DEC 4, 2006

The exudate is further reduced. PolyMem Silver dressings are now changed every three days without additional wound cleansing. Still pain free.




DEC 20, 2006

The wound continues to fill in rapidly. It is staying warm under the PolyMem dressing. Dressing changes are quick and easy: no rinsing.



JAN 2, 2007

The wound is almost closed. The patient has had no antibiotics and no pain medications. She has had no infection or pain since the 2nd week of the new treatment.



JAN 15, 2007

The wound is closed. For the next month, PolyMem dressings, changed every 5 – 6 days, help prevent pain and protect the fragile new soft tissue.

RESULTS

The cleansing provided by the PolyMem dressings supported rapid autolytic debridement and completely eliminated the need for wound cleansing at dressing changes, helping maintain an appropriate wound temperature and promoting quick healing. Within two weeks the woman's pain was completely gone without the need for analgesics. At four weeks wound cultures were negative. At ten weeks the wound was shallow, so the PolyMem Wic Silver cavity filler under PolyMem Silver dressing was discontinued. The wound closed completely in only four months.

CONCLUSION

Healing was remarkable considering the patient's debilitated state. The many unique properties of PolyMem dressings directly influenced the rapidity of the healing and freedom from pain this elderly diabetic woman experienced.

OBJECTIVES

1. Discuss the advantages of PolyMem dressings having a built-in wound cleanser, which minimizes disruption to the wound bed tissue and maintains desirable wound temperature.
2. Note the quick healing with PolyMem dressings, which dramatically decreases costs and provides comfort and convenience to the patient.
3. Consider the role that PolyMem dressing's support of autolytic debridement played in healing this wound.

This case study was unsponsored. Ferris Mfg. Corp. contributed to this poster presentation.



Ferris Mfg. Corp.

5133 Northeast Parkway | Fort Worth, TX 76106 USA

Toll Free USA: 1.800.POLYMEM (765.9636) | International: +1 817.900.1301

Web site: www.PolyMem.com

BIBLIOGRAPHY:

1. Beitz AJ, Newman A, Kahn AR, Ruggles T, Eikmeier L. A polymeric membrane dressing with antinociceptive properties: analysis with a rodent model of stab wound secondary hyperalgesia. *J Pain*. 2004 Feb;5(1):38-47.
2. Hess CT. *Clinical Guide: Skin and Wound Care* 6th Ed. Lippincott Williams & Wilkins. Ambler, PA. 2008;228-231,234-5,343-350.
3. Foresman PA, Etheridge CA, Rodeheaver G. A wound dressing evaluation of partial-thickness rat wounds. *SAWC Health Management Pub*. 1991; Annual Meeting Power Point Presentation.
4. Yastrub DJ. Relationship between type of treatment and degree of wound healing among institutionalized geriatric patients with stage II pressure ulcers. *Care Manag J*. 2004 Winter;5(4):213-8.
5. Burd A, Kwok CH, Hung SC, Chan HS, Gu H, Lam WK, Huang L. A comparative study of the cytotoxicity of silver-based dressings in monolayer cell, tissue explant, and animal models. *Wound Repair and Regeneration* 2007;15:94-104.
6. Clay CS, Chen, WY. Wound pain: the need for a more understanding approach. *Journal of Wound Care*. 2005;14:4,181-184.
7. Fleck, CA. Managing wound pain: today and in the future. *Advances in Skin and Wound Care* 2007;20:3,138-145.
8. Kim Y, Lee S, Hong S, Lee H, Kim E. The effects of polymem on the wound healing. *J Korean Soc Plast Reconstr Surg* 1999;109:1165-1172.

ORIGINAL POSTER PRESENTED AT*:

- 17th Conference of the European Wound Management Association. Poster#135. May 2 - 4, 2997. Glasgow, Scotland.
- 3rd Congress of the World Union of Wound Healing Societies. Poster #PF449. June 4 - 8, 2008. Toronto, Ontario Canada.
- NPUAP 11th Biennial Conference. Poster #40. February 27 - March 1, 2000. Arlington, VA USA.

* This version has been modified from the original; it reflects PolyMem branding.

PolyMem, PolyMem Silver, PolyMem Wic, PolyMem Wic Silver, PolyMem Wic Silver Rope, PolyMem Max, Max, PolyMem Max Silver, Shapes, Shapes by PolyMem, The Shape of Healing, The Pink Dressing, SportsWrap, SportsWrapST, More Healing • Less Pain, interlocking circles design, PolyMem For Sports, Not too Loose... Not too Tight... Just Right!, Ferris and FMC Ferris and design are marks owned by or licensed to Ferris. The marks may be registered or pending in the US Patent and Trademark Office and in other countries. Other marks are the property of their respective owners.