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



EN

PolyMem® Shapes® Tube Dressing

1. The dressing is shaped to snugly fit a tracheostomy or gastrostomy tube.
2. Place the dressing with printed side away from the skin.
3. Surround the tube with the dressing and lay flat against the skin.
4. Position the dressing around the gastrostomy tube in a similar matter.
5. Add medical fixation tape to secure dressing if required.

DESCRIPTION:

Shapes by PolyMem® Tube was created to provide a dressing shaped to accommodate a tube site.

Shapes by PolyMem® Wound Dressings are made of moisture-loving polyurethane matrix that contains (1) safe wound cleanser, (2) glycerin moisturizer, and (3) superabsorbents. Dressings have a thin polyurethane backing which breathes and helps keep out contamination. Dressings may be used as a primary, secondary, or combined primary and secondary dressing.

Drug-free PolyMem dressings help to:

- effectively manage and heal wounds
- absorb fluid and provide a moist healing environment
- minimize and relieve persistent and procedural wound pain
- reduce edema, bruising, and the spread of inflammation into surrounding undamaged tissues

INDICATIONS FOR ALL HUMAN BEINGS:

Under the direction of a healthcare professional, for the management of:

- Tracheostomy tube sites
- Drainage tube sites
- Gastrointestinal tube sites
- Feeding tube sites

PRECAUTIONS:

- For external use only
- PolyMem dressings are not compatible with oxidizing agents (hydrogen peroxide or hypochlorite solutions)
- Prior to performing imaging procedures, determine if the dressing should be removed
- Avoid contact with electrodes or conductive gels

- Topical treatments are not recommended in conjunction with PolyMem dressings
- Be alert for signs of infection or maceration
- Do not use and discontinue use on people who show signs of sensitivity, irritation, or allergy from the dressing or its materials
- Intended for one use on one person. Do not reuse. Reuse may degrade properties, spread contamination, and increase risk of infection.

PolyMem dressings may be used when signs of infection are present if proper medical treatment that addresses the cause of the infection has been implemented.

Store at room temperature ≤30°C (86°F). Excursions permitted ≤55°C (131°F) but should be minimized.

Disposal – as is typical in the facility and according to local requirements:

- used devices - biohazardous waste
 - unused or expired devices - non-hazardous waste
- Dimensions** are approximate.

APPLICATION

Initial Dressing Application:

1. Prepare the tube site according to protocol or as directed by a physician or other ordering clinician.
2. Apply the dressing film side out and/or printed side out. If desired, the dressing can be cut to shape prior to application. Gently slide the Tube dressing so that it fits neatly and easily around the tube.
3. Medical fixation tape may also be applied to secure the dressing to the surrounding skin.
4. Topical treatments are not recommended in conjunction with Shapes® dressings.
5. Keep the dressing DRY and in place when bathing. Change the dressing if it gets wet.

Dressing Change

1. Remove the dressing when wound fluid or secretions, absorbed into the dressing, are in contact with intact skin, or when good practice dictates.
2. Cleanse the tube area, if necessary.
3. Apply a new dressing.

Excessive cleaning may injure regenerating tissue and delay wound healing.


Observations during use or change:

Be alert for signs and symptoms, such as, increased redness, inflammation, pain, warmth, odor, whitening, fragility, excessive tissue formation, or lack of healing. Consult a professional for problem situations, such as, infections, irritations, sensitivity, allergy, maceration, or hypergranulation.

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Ferris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.

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© 2017, 2019 Ferris Mfg. Corp.® or ™ indicate trademark