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
EN


## PolyMem® Finger/Toe Dressing

*The dressing provides a soothing wrap shaped to fit.*

**Size Selection:**  
Use the image on the outside of the box or pouch to help choose the correct size of the dressing (1, 2, 3, 4, or 5).

**Application:** The dressing was designed for easy application; simply roll onto the finger.

### The grid and words should be on the outside of the dressing.

The dressing should be snug, but maintain comfort. If it feels too tight, remove and apply a larger size. Remove by rolling off, opposite of application. The dressing may be applied in different ways. For toes the dressing may be cut along the sides creating flaps that can be laid upon the top and bottom of the foot. These can be secured with tape, or the dressing may be cut to form a ring or sleeve over the injured portion of the finger.

### DESCRIPTION:

PolyMem polymeric membrane dressings are made of moisture-loving polyurethane matrix that contains (1) safe wound cleanser, (2) glycerol moisturizer, and (3) superabsorbents. Dressings have a thin polyurethane backing which breathes and keeps out contamination. Dressings may be used as primary, secondary, or combined primary and secondary dressing.

Drug-free PolyMem dressings help to:

- effectively manage and heal wounds
- absorb fluid and provide a moist healing environment
- minimize and relieve persistent and procedural wound pain
- reduce edema, bruising and the spread of inflammation into surrounding undamaged tissues

### INDICATIONS FOR ALL HUMAN BEINGS:

**Over-the-counter**, for management of minor wounds, such as:

- Abrasions • Lacerations • Cuts • Burns
- **Under the direction of a healthcare professional**, for the management of serious wounds, such as most partial and full thickness wounds, including:
  - Pressure ulcers (stages I-IV) • Acute wounds • Skin tears • 1st and 2nd degree burns • Leg ulcers • Diabetic ulcers • Surgical wounds
  - Venous stasis ulcers • Donor and graft sites • Dermatologic disorders • Vascular ulcers

### PRECAUTIONS:

- For external use only • PolyMem dressings are not compatible with oxidizing agents (hydrogen peroxide or hypochlorite solutions) • Prior to performing imaging procedures, determine if the dressing should be removed • Avoid contact with electrodes or conductive gels
- Topical treatments are not recommended in conjunction with

PolyMem dressings • Be alert for signs of infection or maceration • Do not use and discontinue use on people who show signs of sensitivity, irritation, or allergy from the dressing or its materials • Intended for one use on one person. Do not reuse. Reuse may degrade properties, spread contamination, and increase risk of infection

**PolyMem dressings may be used when signs of infection are present if proper medical treatment that addresses the cause of the infection has been implemented.**

**Store** at room temperature ≤30°C (86°F). Excursions permitted ≤55°C (131°F) but should be minimized.

**Disposal** – as is typical in the facility and according to local requirements:

- used devices - biohazardous waste
  - unused or expired devices - non-hazardous waste
- Dimensions** are approximate.

### DIRECTIONS:

#### Initial Application:

1. Prepare the injured area according to protocol or as directed by a physician or other ordering clinician. Clean the skin of dirt, topical medications and lotions. Residual moisture on the skin is beneficial.
2. Apply the dressing as outlined. A PolyMem dressing should also cover any inflamed or damaged area surrounding the injury. Dressings may be cut to accommodate the injury or body part.
3. Apply film side out (so printing is visible).

#### Helpful Hints:

- Secure dressing with water-proof roll tape to reduce the chance of the inside of the dressing becoming contaminated.
- Outline the approximate wound margin on the outside of the dressing to help determine when to change.

#### Use:

1. A dramatic increase in wound fluid may be observed during the first few days. This is not uncommon, and indicates that the dressing is working.
2. Keep the dressing dry and in place when bathing or performing hand-hygiene. Change the dressing if it gets wet.

#### Dressing Change:

1. For an exuding wound, fluid will become visible through the top. Change before fluid reaches the wound margin, or when good practice dictates, or after no more than 7 days. Change immediately if fluid reaches the edge of the membrane pad.

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2. More frequent changes may be needed due to a compromised immune system, diabetes, infection, hypergranulation, macerated tissue, or when desiring to facilitate the removal of non-viable tissue.
3. Gently remove the dressing. The PolyMem pad will generally not adhere to the wound bed, usually assuring pain-free dressing changes. Inspect the dressing and wound to assure removal of all material.
4. If coagulated blood or dried exudate has caused the dressing to adhere to skin or hair, use a small amount of saline or water to slowly gently soften and loosen. Clean the intact skin around the wound.
5. Manage dressing changes according to protocol or as directed by a physician or other ordering clinician. In most cases, when using PolyMem, there is no need to disturb or cleanse a wound during changes unless the wound is infected or contaminated. PolyMem is designed to continuously clean the wound and does not leave residue that needs to be removed.
- Excessive cleaning may injure regenerating tissue and delay wound healing.**
6. Apply a new dressing.

#### Observations during use or change:

- During the normal healing process, the wound may appear to become larger during the first few dressing changes because the body is naturally removing non-viable tissue. If the wound continues to become larger, consult an expert.
- PolyMem dressings help support new blood vessel formation so it is not uncommon to see blood-stained wound fluid and dressings during changes.
- Be alert for signs and symptoms, such as, increased redness, inflammation, pain, warmth, odor, whitening, fragility, excessive tissue formation, or lack of healing. Consult a professional for problem situations, such as infection, irritation, sensitivity, allergy, maceration, or hypergranulation.

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Ferris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.

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© 2025 Ferris Mfg. Corp.® or ™ indicate trademark