
<p>Do not use if pouch has been damaged or opened. Nicht verwenden, wenn der Beutel beschädigt oder geöffnet ist. Använd inte skadad eller öppnad förpackning. Ikke bruk hvis forpakningen er skadet eller åpenet. Non utilizzare se la busta è stata aperta o danneggiata. Ne pas utiliser si l'emballage du pansement a été ouvert ou est endommagé. No utilizar si el estuche está dañado o abierto. Eğer poş zarar görmüş veya açılmış ise kullanmayın.</p>	<p>Keep away from sunlight and heat. Vor Sonnenlicht und Hitze schützen. Håll borta från solljus och värme. Ikke utsett produktet for sollys og varme. Tenere lontano dalla luce solare e dal calore. Garder à l'abri des rayons du soleil et de la chaleur. Mantener fuera del alcance de la luz solar y el calor. Güneş ışığı ve ısıdan koruyun.</p>
<p>Store ≤30°C (86°F). Bei ≤30°C (86°F) lagern. Förvara i ≤30°C (86°F). Lagre ved ≤30°C (86°F). Conservare a ≤30°C (86°F). Conserver à une température ≤30°C (86°F). Consevar a ≤ 30°C (86°F). ≤30°C (86°F) ve altındaki sıcaklıklarda saklayın.</p>	<p>Not made with natural rubber latex. Nicht aus Naturkautschuklatex hergestellt. Produkten är inte tillverkad av naturgummlatex. Ikke laget av naturlig lateks. Privo di lattice di gomma naturale. Ne contient pas de latex naturel. No están hechos de látex de caucho natural. Doğal kauçuk lateksten üretilmemiştir.</p>
	<p>Single use. Do not re-use. Zum Einmalgebrauch. Nicht wiederverwenden. För engångsbruk. Återanvänd inte. Til engangsbruk. Ikke bruk om igjen. Monouso. Non riutilizzare. À usage unique. Ne pas réutiliser. Desechables. No reutilizar. Tek kullanımlıdır. Tekrar kullanmayın.</p>

by **PolyMem**[®]
SPORTSWRAP[®]
ST

REF 9120

2.5" x 12'
(6.3 cm x 3.6 m)

FERRIS MFG. CORP.
5133 Northeast Parkway
Fort Worth, TX 76106-1822 U.S.A.
Toll Free U.S.A.: 800-765-9636
International: +1 817-900-1301
www.polymem.com
© 2026 Ferris Mfg. Corp.® or ™ indicate trademark

Emergo Asia Pacific Pty Ltd
201 Sussex Street
Darling Park, Tower II, Level 20
Sydney, NSW 2000 Australia

EU REP

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

AW-0573_R02 0326
FMC-5910_R7 0326

EN

ROLL

PRODUCT DESCRIPTION
SportsWrap ST provides the same benefits as SportsWrap, but in a slimmer, thinner profile. It may be:

- applied under tape, wraps, or splints
- used during rest, rehab and recovery; and during activity, training and competition SportsWrap ST roll is a patented hydrophilic polyurethane membrane material.

INDICATIONS
For the local management of swelling, inflammation, bruising, stiffness or pain due to acute or chronic injury when the skin is intact.

PRECAUTIONS
Consult your physician for all serious injuries, and if symptoms persist after wearing or increase during use. Not sterile. Not for broken skin or open wounds. For external use only.

PREPARATION
Clean the skin of dirt, topical medications and lotions. Residual moisture on skin is OK.

APPLICATION
Apply as soon as possible after injury. For chronic injury management apply prior to physical activity.
1. Unroll the green membrane wrap in snug contact with skin.
2. Wrap in a continuous motion overlapping only slightly. Cut or tear when complete.
3. Secure with athletic or adhesive tape, elastic wrap or garment. Wrap firmly for close contact with the skin but apply gently to maintain comfort. The green membrane must maintain intimate contact with the skin to work effectively. Return the unused portion to the pouch.

WEARING
1. Avoid wetting. Remove before bathing or showering.
2. Wear during exercise. After exercise, discard and replace.
3. Saturated material cannot be dried out and reused.

STORAGE
Keep dry and store at room temperature and humidity.

SportsWrap ST rolls compared to standard SportsWrap rolls and wraps:

Do both products provide the same great performance?
Both dressings provide the same key benefits to the athlete and the athletic trainer:
• help to manage strains, sprains, bruises and contusions where the skin is intact
• help to reduce pain, swelling, bruising and inflammation
• protect the injury
• adjunct to and synergistic with most physical modalities and other therapies
• simple and convenient to use
• use 24/7 after injury
• wear during rest and sleep
• long-lasting, drug-free pain relief
• get the athlete back in the game faster!

This combination of benefits is available only with ST and SportsWrap products!

What is different about ST?
The material is slim and thin for ease of use.

What experience have athletic trainers had with ST?
ST rolls have been used successfully by athletic trainers of elite professional, international and college teams since ST's introduction in 2010.

When to use ST?
Use especially for short-term response (after injury, or during initial recovery), or to help manage the pain and swelling that accompanies chronic injuries.

How to use ST after injury?
1. Apply ST similar to other wrapping materials. Also protect the injury with tape, splint or wrap.
2. On the following day, remove and check the injury. Discard. Apply again.
3. Repeat this cycle several times.

How to use ST for return to practice?
1. May be used similar to tape underwrap to provide the benefits of SportsWrap during activity.
2. After practice, remove and check the injury and athlete. Discard and apply again.
3. May be combined with SportsWrap to continue benefits while athlete is at rest.

When to use standard SportsWrap?
1. Use SportsWrap for immediate management of acute sprains, strains and contusions; and to assist with pain and inflammation control of chronic injuries.
2. Applying SportsWrap rolls and wraps has the ease, adjustability and convenience of a conventional elastic bandage. Therefore, the athlete may be instructed in use, giving the athlete options for when and how to wear the SportsWrap products, resulting in better compliance with the overall recovery program.
3. The same SportsWrap roll or wrap may be used for multiple days. The athlete should place the SportsWrap in a plastic bag for protection and convenience when not being worn.

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Ferris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.

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