

PolyMem® Ordering Information


PolyMem and PolyMem Silver Non-Adhesive Dressings

PolyMem Non-Adhesive				
REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5022	1.8" x 1.8" (4.7cm x 4.7cm)	100 per case	5 bx / 20 ea	A6209
5033	3" x 3" (8cm x 8cm)	60 per case	4 bx / 15 ea	A6209
5044	4" x 4" (10cm x 10cm)	60 per case	4 bx / 15 ea	A6209
5055	5" x 5" (13cm x 13cm)	30 per case	2 bx / 15 ea	A6210
5077	6.5" x 7.5" (17cm x 19cm)	30 per case	2 bx / 15 ea	A6210
5124	4" x 12.5" (10cm x 32cm)	12 per case	1 bx / 12 ea	A6210
5244	4" x 24" (10cm x 61cm)	8 per case	2 bx / 4 ea	A6211
5824	8" x 24" (20cm x 60cm)	2 per case	1 bx / 2 ea	A6211
PolyMem Silver Non-Adhesive				
1022	1.8" x 1.8" (4.7cm x 4.7cm)	100 per case	5 bx / 20 ea	A6209
1044	4.25" x 4.25" (10.8cm x 10.8cm)	30 per case	2 bx / 15ea	A6210
1077	6.5" x 7.5" (17cm x 19cm)	15 per case	1 bx / 15ea	A6211
1124	4.25" x 12.5" (10.8cm x 32cm)	12 per case	1 bx / 12ea	A6211

PolyMem®

PolyMem® SILVER

The optimal choice when cutting dressings to size and taping is your preferred method of application. Use as a combined primary and secondary dressing or as a secondary dressing.




PolyMem Film Adhesive Dressings

PolyMem Film Adhesive				
REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
203	1" x 1" (2.5cm x 2.5cm) Pad 2" x 2" (5cm x 5cm) Adhesive	100 per case	5 bx / 20 ea	A6212
405	2" x 3" (5cm x 8cm) Pad 4" x 5" (10cm x 13cm) Adhesive	60 per case	4 bx / 15 ea	A6212
606	3.5" x 3.5" (9cm x 9cm) Pad 6" x 6" (15cm x 15cm) Adhesive	60 per case	4 bx / 15 ea	A6212
3412	2" x 10" (5cm x 26cm) Pad 4" x 12.5" (10cm x 32cm) Adhesive	12 per case	1 bx / 12 ea	A6213

PolyMem®

The ultra-thin, adhesive-coated urethane film, together with the PolyMem membrane, will provide a breathable yet tough, protective, water-resistant dressing of unsurpassed versatility.




PolyMem Max and PolyMem Max Silver Dressings

PolyMem Max				
REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5035	3" x 3" (7.6cm x 7.6cm)	40 per case	4 bx / 10 ea	A6209
5045	4.5" x 4.5" (11cm x 11cm)	20 per case	2 bx / 10 ea	A6210
5088	8" x 8" (20cm x 20cm)	10 per case	2 bx / 5 ea	A6211
PolyMem Max Silver				
1045	4" x 4" (10cm x 10cm)	16 per case	2 bx / 8 ea	A6209
1088	8" x 8" (20cm x 20cm)	10 per case	2 bx / 5 ea	A6211

PolyMem® MAX

This configuration is a thicker PolyMem formulation for wounds with heavier drainage or when longer wear time is desired. Use as a combined primary and secondary dressing or as a secondary dressing.




PolyMem Wic, PolyMem Wic Silver & PolyMem Wic Silver Rope Wound Fillers

PolyMem Wic				
REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5733	3" x 3" (8cm x 8cm) 4 grams	40 per case	4 bx / 10 ea	A6215
5712	3" x 12" (8cm x 30cm) 16 grams	12 per case	1 bx / 12 ea	A6215
5788	8" x 8" (20cm x 20cm) 28 grams	10 per case	2 bx / 5 ea	A6215
PolyMem Wic Silver				
1331	1" x 3" (2.5cm x 7.6cm) 1.3 grams	70 per case	5 bx / 14 ea	A6215
1333	3" x 3" (8cm x 8cm) 4 grams	20 per case	2 bx / 10 ea	A6215
1788	8" x 8" (20cm x 20cm) 28 grams	10 per case	2 bx / 5 ea	A6215
PolyMem Wic Silver Rope				
1814	0.4" x 14" (1cm x 35cm) 3 grams	12 per case	2 bx / 6 ea	A6215

PolyMem® WIC

With no film backing, Wic expands within the wound cavity to fill dead space, keeping the wound bed clean and moisturized while absorbing excess fluid. Wic Rope has a reinforced mesh to ensure one-piece removal.




PolyMem Cloth Adhesive Dressings

PolyMem Cloth Adhesive				
REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
7203	1" x 1" (2.5cm x 2.5cm) Pad 2" x 2" (5cm x 5cm) Adhesive	100 per case	5 bx / 20 ea	A6212
7405	2" x 3" (5cm x 8cm) Pad 4" x 5" (10cm x 13cm) Adhesive	60 per case	4 bx / 15 ea	A6212
7606	3.5" x 3.5" (9cm x 9cm) Pad 6" x 6" (15cm x 15cm) Adhesive	60 per case	4 bx / 15 ea	A6212
7031	1" x 1" (2.5cm x 2.5cm) Pad 1" x 3" (2.5cm x 8cm) Adhesive	100 per case	5 bx / 20 ea	A6212
7042	2" x 1.5" (5cm x 4cm) Pad 2" x 4" (5cm x 10cm) Adhesive	100 per case	5 bx / 20 ea	A6212
PolyMem Silver Cloth Adhesive				
1766	3.5" x 3.5" (9cm x 9cm) Pad 6" x 6" (15cm x 15cm) Adhesive	30 per case	2bx / 15 ea	A6212

PolyMem®


A breathable, adhesive-coated cloth backing that provides exceptional comfort for your patient.



PolyMem Finger/Toe Dressings

PolyMem	Silver	Finger / Toe			
REF	REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
4401	1401	#1 small; circumference 1.8" - 2.2" (46.7mm - 57.0mm)	30 per case	5 bx / 6 ea	A6209
4402	1402	#2 medium; circumference 2.2" - 2.6" (57.0mm - 67.2mm)	30 per case	5 bx / 6 ea	A6209
4403	1403	#3 large; circumference 2.6" - 3.0" (67.2mm - 77.4mm)	30 per case	5 bx / 6 ea	A6209
4404	1404	#4 extra large; circumference 3.0" - 3.4" (77.4mm - 87.6mm)	30 per case	5 bx / 6 ea	A6209
4405	1405	#5 XXL; circumference 3.4" - 3.8" (87.6mm - 97.8mm)	30 per case	5 bx / 6ea	A6210

This configuration is specially designed for use on fingers and toes. Appropriate for use as a combined primary and secondary dressing or as a secondary dressing.




Shapes by PolyMem Oval & Oval Silver Dressings

PolyMem	Silver	Shapes by PolyMem Oval			
REF	REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
8023	1823	1" x 2" (2.5cm x 5cm) Pad 2" x 3" (5cm x 7.6cm) Adhesive	100 per case	5 bx / 20 ea	A6212
8053	1853	3" x 2" (7.6cm x 5cm) Pad 5" x 3.5" (12.7cm x 8.8cm) Adhesive	30 per case	2 bx / 15 ea	A6212
8086	1886	4" x 5.7" (10.1cm x 14.6cm) Pad 6.5" x 8.2" (16.5cm x 20.9 cm) Adhesive	20 per case	2 bx / 10 ea	A6213

Shapes by PolyMem®

A variety of oval sizes offer a secure fit for wounds on elbows, knees and other contoured areas of the body where the dressing needs to bend and conform.



Shapes by PolyMem Sacral Dressings


PolyMem	Silver	Shapes by PolyMem Sacral			
REF	REF	Dimensions	PER CASE	PACKAGING	HCPCS Code
3709	1709	4.5" x 4.7" (11.4cm x 12.0cm) Pad 7.2" x 7.8" (18.4cm x 20cm) Adhesive	20 per case	2 bx / 10 ea	A6212



The heart-shaped dressing is designed specifically for use in the sacral region.

Shapes by PolyMem Tube Dressings

Shapes by PolyMem Tube				
REF	PAD DIMENSIONS	PER CASE	PACKAGING	HCPCS Code
5333	2.75" x 2.75" (7cm x 7 cm)	30 per case	2 bx / 15 ea	A6209
5335	3.5" x 3.5" (9cm x 9 cm)	30 per case	2 bx / 15 ea	A6209



This configuration is pre-cut to enhance patient comfort and help clinicians easily apply the dressing to tube sites. Unlike other dressings used on tube sites, PolyMem dressings do not shed fibers.



Ferris Mfg. Corp. and its representatives provide general information as a courtesy and for the awareness of customers. Ferris does not guarantee coverage or payment for products or accuracy of HCPCS codes associated with REF numbers (thought to be current 1-1-2016.) Ferris information is not intended to increase or maximize reimbursement by any payer, and does not recommend how products should be billed to Medicare or other health care insurers. The Provider bears responsibility for coding, satisfying requirements, and obtaining reimbursement. The Provider's records and billing should accurately describe the product and amount that was clinically necessary to treat the condition and the amount provided to the patient. Coverage, coding, and payment policies and guidelines are complex, can vary from one carrier or region to another, and change. Providers should consult with counsel, a reimbursement specialist, and CMS web pages. Providers should check with their local carriers or intermediaries often and consider the DME PDAC DMECS Application and calling the local DME MAC.

This document is meant for general informational purposes only. See individual product literature for specific indication and instructions for use.

PolyMem, PolyMem Silver, PolyMem Wic, Wic, PolyMem Wic Silver, PolyMem Wic Silver Rope, PolyMem Max, Max, PolyMem Max Silver, Shapes, Shapes by PolyMem, The Shape of Healing, The Pink Dressing, SportsWrap, SportsWrapST, More Healing • Less Pain, interlocking circles design, PolyMem For Sports, Not too Loose... Not too Tight... Just Right!, Ferris and FMC Ferris and design are marks owned by or licensed to Ferris. The marks may be registered or pending in the US Patent and Trademark Office and in other countries. Other marks are the property of their respective owners.

© 2016 Ferris Mfg. Corp.

MKL-698 REV0 0116



MADE IN USA